

Your partner for
improved hearing.

powerone)))



www.powerone-batteries.com

power one –
proven reliability,
proven performance.

Your requirements

Do you need hearing aid batteries that function safely and reliably, that you can rely on in any situation and that offer extraordinary performance? Then try **power one** and experience the benefits for yourself.

power one manufactures its products in the world's largest and most modern production facility for hearing aid batteries in Ellwangen, Germany.



- **power one** is a premium battery with a longer life – even for digital and high-power hearing aids.
- **power one** is reliable and safe – even in extreme environmental conditions such as high humidity.
- **power one** lets you to take full advantage of your hearing aid's capacity at all times, ensuring brilliant sound quality.
- **power one** is made in Germany and is subject to the strictest quality controls.
- **power one** offers practical, convenient and easy-to-handle packaging.

power one –
simply comfortable,
simply practical.

Your expectations

Do you need packaging that is easy to handle and keeps your batteries safe?

power one has turned your hope into reality:

- Practical rotating tray lets you remove the batteries easily and conveniently
- Special quality seal ensures maximum product safety
- Safe storage using recloseable clip
- Easy-to-see colour coding lets you identify your battery type immediately
- Environmentally friendly PET and cardboard materials



safety comfort pack™

The practical rotating tray lets you easily remove and re-insert the cells, conveniently and securely. Two layers of cardboard ensure that no batteries fall out.

comfort tab

The extra-long tab lets you conveniently insert the battery into the hearing aid.

power seal

The quality seal guarantees that the packaging has not been opened.



Better performance
through proper handling.

Activation

power one is only activated once you remove the comfort tab. Therefore, you should wait until just before inserting the battery into the hearing aid before pulling off the tab, as the battery uses oxygen in the air to operate.

Storage

Store your battery in a dry place. The ideal storage temperature is between 10°C and 30°C.



Disposal

Never leave an empty hearing aid battery in the hearing aid longer than necessary. Dispose of the battery in an environmentally sound manner or return the used battery to your hearing aid specialist.

Protection

Never allow the battery to come into contact with metal objects or discharge will occur.

Store your battery in its intended packaging.

www.powerone-batteries.com